

<b>Shri Manibhai Virani &amp; Smt. Navalben Virani Science College</b> Yogidham Gurukul, Kalawad Road, Rajkot – 360005, Gujarat, India		
Doc. no. <b>VSC/CIF/SOP/001</b>	Doc. name: <b>Standard Operating Procedure for GCMS</b>	Location : <b>CIF Lab-111</b>
Model : <b>GC(7890B), MS(5977A)</b> Sr.No : <b>CN14153010</b>	Make: <b>Agilent Technologies</b>	ID no. <b>VSC/CIF/001</b>
Prepared by: Arjun Borad		Approved by: Dr. K. D. Ladva

## 1. OBJECTIVE

To describe the procedure for operation and maintenance of GCMS

## 2. SCOPE

Applicable to GC-MS system in the Central Instrumentation Facility (CIF) Lab

## 3. REFERENCE & ATTACHMENTS :

- 3.1 User's operating manual supplied by manufacturer.
- 3.2 Attachment – I Log book record.
- 3.3 Attachment –II Maintenance Log

## 4. RESPONSIBILITY

### 4.1 Head of Department

- 4.1.1 To ensure the implementation of this SOP.
- 4.1.2 To initiate repairs and to make alternative arrangements during the pendency of repair.

### 4.2 Application Chemist

- 4.2.1 To ensure proper documentation as per SOP.
- 4.2.2 To maintain the instrument and its accessories.
- 4.2.3 To inform about any failure of instrument or calibration failure to Head and to make a record in maintenance log.

## 5. ABBREVIATIONS

<b>SOP</b>	Standard operating procedure
<b>GCMS</b>	Gas Chromatography Mass Spectroscopy
<b>kPa</b>	Kilopascal
<b>GC</b>	Gas Chromatography

## 6. PROCEDURE

- This section describes following basic operating procedure for GCMS:
  - Gas requirements
  - Instrument Startup & Shutdown

- Maintenance
- An understanding of this topics is crucial for operating GCMS instrument and GCMS Online software.

## 6.1 Gas requirements

- This section describes the carrier gas used by this system. The specifications presented below must be followed to promote safety and maintain instrument performance.
- Carrier gas & Supply Pressure
- Helium (300-980 kPa)
- Hydrogen (Depends on Analytical condition)
- Methane (Depends on Analytical condition)
- Isobutane (Depends on Analytical condition)
- Argon (350-450 kPa)
- Oxygen (below 3 ppm)
- The Pressure and flow rate setting range on GC varies according to supply pressure..
- **Note:** Any air leakage into the gas supply lines will degrade the purity of gas and may affect system performance.

## 6.2 Instrument Startup & Shutdown

- 6.2.1 Turn on the Carrier gas.
- 6.2.2 Turn on the power switches of whole instrument.
- 6.2.3 First start GC, MS respectively.
- 6.2.4 Turn on Computer, Monitor and start window.
- 6.2.5 Double click a GCMS Online.
- 6.2.6 Select the used id and enter the password.
- 6.2.7 Click OK so the program starts up.
- 6.2.8 Close the GCMS Online.
- 6.2.9 Switch OFF the GC, MS respectively.
- 6.2.10 Switch OFF the power switches of whole instrument.

## 6.3 Maintenance

- 6.3.1 Make sure the Column is Connected at both side (Injection port and Detector side)
- 6.3.2 Keep the instrument in a dust, vibrations, spatial noise and corrosive gases free environment.
- 6.3.3 Wipe the instrument to make it dust free before operation.

- 6.3.4 Do not turn off power supply to the main electronic of the instrument and pump at any time except under maintenance.
- 6.3.5 Make sure the continuous supply of carrier gas.
- 6.3.6 Check for leaks, always after performing maintenance.
- 6.3.7 Flush the column thoroughly with appropriate solvent after each sample analysis.
- 6.3.8 Oil in the pump should be monitored to maintain the minimum level.

#### 7. DISTRIBUTION OF SOP:

- Master copy (original): To be kept with Head of Department.
- Reference copy: To be kept with Application Chemist
- Reference copy: For display at instrument



  
Principal  
Shri Manibhai Virani and  
Smt. Navalben Virani Science College  
(Autonomous) Rajkot.